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1. Foley J, Evans D, Blackwell A. Partial caries removal and cariostatic materials in carious primary molar teeth: A randomised controlled clinical trial. Br Dent J. 2004;197(11):697-701.
2. Mertz-Fairhurst EJ, Curtis JW Jr, et al. Ultraconservative and cariostatic sealed restorations: results at year 10. J Am Dent Assoc. 1998;129(1):55-66.
3. Hilton TJ. Keys to clinical success with pulp capping: a review of the literature. Oper Dent. 2009;34(5):615-625.
4. Demarco FF, Tarquinio SB, Jaeger MM, et al. Pulp response and cytotoxicity of 2 dentin bonding agents. Quintessence Int. 2001;32(3): 211-220.
5. Cehreli ZC, Turgut M, Olmez S, et al. Short term human primary pulpal response after direct pulp capping with fourth-generation dentin adhesives. J Clin Pediatr Dent. 2000;25(1):65-71.
6. Kamitakahara M, Kawashita M, Kokubo T, Nakamura T. Effect of polyacrylic acid on the apatite formation of a bioactive ceramic in a simulated body fluid: fundamental examination of the possibility of obtaining bioactive glass-ionomer cements for orthopaedic use. Biomaterials. 2001;22(23):3191-3196.

7. Hebling J, Costa C, Nogueira I, et al. Cytotoxicity of resin-based light-cured liners. Am J Dent. 2009;22(3):137-142.
8. Cannon M, Gerodias I, Vieira EM, et al. Effects of Different Pulp Capping Techniques on Hard Tissue Bridge Formation. European Academy of Pediatric Dentistry. June 2006, Amsterdam, Netherlands.
9. Mineral Trioxide Aggregate, Comprehensive Literature Review, Journal of Endodontics, March 2010.
10. Gandolfi MG, Suh B, Siboni F. Chemical-physical properties of TheraCal LC pulp capping material. Presented at: International Association of Dental Research (IADR). March 18, 2011; San Diego, CA. Abstract #2521.
11. Gandolfi MG, Siboni F, Taddei P, et al. Apatite-forming ability of TheraCal LC pulp capping material. Presented at: International Association of Dental Research (IADR). March, 2011; San Diego, CA. Abstract #2520.
12. Gandolfi MG, Siboni F, Taddei P, et al. Chemical-physical properties of TheraCal LC, a novel light-curable MTA-like material for pulp capping. International Endodontic Journal. 2012.

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Chemical–physical properties of TheraCal, a novel light-curable MTA-like material for pulp capping

M. G. Gandolfi, F. Siboni & C. Prati

Laboratory of Biomaterials and Oral Pathology, Department of Odontostomatological Sciences, University of Bologna, Bologna, Italy

Received 14 July 2011; accepted 29 December 2011. International Endodontic Journal

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- Bogen G, Kim JS, Bakland LK (2008) Direct pulp capping with mineral trioxide aggregate. An Observational Study. *Journal of American Dental Association* 139, 305–15.
- Clapham DE (1995) Calcium signaling. *Cell* 80, 259–68.
- Daneš G, Dammaschke T, Gerth HUV, Zandbiglari T, Schafer E (2006) A comparative study of selected properties of ProRoot mineral trioxide aggregate and two Portland cements. *International Endodontic Journal* 39, 213–9.
- Desai S, Chandler N (2009) Calcium hydroxide-based root canal sealers: a review. *Journal of Endodontics* 39, 415–22.
- Devito LK, Ortega AI, Haiter-Neto F (2006) Effect of the storage in water on the radiopacity of calcium hydroxide cements. *Brazilian Journal of Oral Sciences* 5, 958–62.
- Dougherty EW (1962) Dental cement material. US patent 3047408.
- Estrela C, Holland R (2003) Calcium hydroxide study based on scientific evidences. *Journal of Applied Oral Sciences* 11, 269–82.
- European Society of Endodontology (2006) Quality guidelines for endodontic treatment: consensus report of the European Society of Endodontology. *International Endodontic Journal* 39, 321–30.
- Furey A, Hjeltnaug J, Lobner D (2010) Toxicity of Flow Line, Durafill VS, and Dycal to dental pulp cells: effects of growth factors. *Journal of Endodontics* 36, 1149–53.
- Gandolfi MG, Pagani S, Perut F, Ciapetti G, Baldini N, Prati C (2008) Innovative silicate-based cements for endodontics: a study of osteoblast-like cell response. *Journal of Biomedical Materials Research* 86(A), 477–86.
- Gandolfi MG, Ciapetti G, Perut F et al. (2009) Biomimetic calcium-silicate cements aged in simulated body solutions. Osteoblasts response and analyses of apatite coating. *Journal of Applied Biomaterials and Biomechanics* 7, 160–70.
- Gandolfi MG, Taddei P, Tinti A, Prati C (2010a) Apatite-forming ability of ProRoot MTA. *International Endodontic Journal* 43, 917–29.
- Gandolfi MG, Van Landuyt K, Taddei P, Modena E, Van Meerbeek B, Prati C (2010b) ESEM-EDX and Raman techniques to study ProRoot MTA and calcium-silicate cements in wet conditions and in real-time. *Journal of Endodontics* 36, 851–7.
- Gandolfi MG, Taddei P, Tinti A, Dorigo De Stefano E, Rossi PL, Prati C (2010c) Kinetics of apatite formation on a calcium-silicate cement for root-end filling during ageing in physiological-like phosphate solutions. *Clinical Oral Investigations* 14, 659–68.
- Gandolfi MG, Ciapetti G, Taddei P et al. (2010d) Effect of ageing on bioactivity and in vitro biological properties of calcium-silicate cements for dentistry. *Dental Materials* 26, 974–92.
- Gandolfi MG, Shah SN, Feng R, Prati C, Akintoye SO (2011a) Biomimetic calcium-silicate cements support differentiation of human orofacial bone marrow stromal cells. *Journal of Endodontics* 37, 1102–8.
- Gandolfi MG, Taddei P, Siboni F, Modena E, Ginebra MP, Prati C (2011b) Fluoride-containing nanoporous calcium-silicate MTA cements for endodontics and oral surgery: early fluorapatite formation in a phosphate-containing solution. *International Endodontic Journal* 44, 938–49.
- Gandolfi MG, Taddei P, Tinti A, De Stefano Dorigo E, Prati C (2011c) Alpha-TCP improves the apatite formation ability of calcium-silicate hydraulic cement soaked in phosphate solutions. *Materials Science Engineering C* 31, 1412–22.
- Gandolfi MG, Taddei P, Siboni F, Modena E, Ciapetti G, Prati C (2011d) Development of the foremost light-curable calcium-silicate MTA cements root-end in oral surgery. Chemical-physical properties, bioactivity and biological behaviour. *Dental Materials* 27, e134–57.
- Gandolfi MG, Taddei P, Siboni F, Modena E, Prati C (2011e) Biomimetic remineralization of human dentine using promising innovative calcium-silicates hybrid "smart" materials. *Dental Materials* 27, 1055–69.
- Golberg M, Smith AJ (2004) Cells and extracellular matrices of dentin and pulp/a biological basis for repair and tissue engineering. *Critical Reviews in Oral Biology & Medicine* 15, 13–27.
- Hlebling J, Lessa FCR, Nogueira I, Carvalho RM, de Souza Costa CA (2009) Cytotoxicity of resin-based light-cures linal cements. *Journal of Dental Research* 87B, 470.
- Horsted-Bindslev P, Lovshall H (2002) Treatment outcome of vital pulp treatment. *Endodontic Topics* 2, 24–34.
- Hosoya N, Takahashi G, Arai T, Nakamura J (2001) Calcium concentration and pH of the periapical environment after applying calcium hydroxide into root canals in vitro. *Journal of Endodontics* 27, 343–6.
- Islam I, Chng HK, Yap AUJ (2006) Comparison of the physical and mechanical properties of MTA and Portland cement. *Journal of Endodontics* 32, 193–7.
- ISO 4049 (2010) Dentistry-polymer-based filling, restorative and luting materials. Geneva: International Organization for Standardization ISO 4049.
- ISO 6876 (2002) Dental root canal sealing materials. Geneva: International Organization for Standardization ISO 6876.
- Jung GY, Park YJ, Han JS (2010) Effects of HA released calcium ion on osteoblast differentiation. *Journal of Materials Science: Materials in Medicine* 21, 1649–54.
- Lazic S (1995) Microcrystalline hydroxyapatite formation from alkaline solutions. *Journal Crystal Growth* 147, 147–54.
- Lopez-Cazaux S, Blutescu G, Mague D, Lieubeau B, Guicheux J, Alliot-Licht B (2006) Culture medium modulates the behaviour of human dental pulp-derived cells: technical note. *European Cells & Materials Journal* 11, 35–42.
- Maeno S, Niki Y, Matsumoto H et al. (2005) The effect of calcium ion concentration on osteoblast viability, proliferation and differentiation in monolayer and 3D culture. *Biomaterials*, 26, 4847–55.32.Meyer JL, Eanes ED (1978) A thermodynamic analysis of the amorphous to crystalline phosphate transformation. *Calcified Tissue Research* 25, 59–68.
- Min KS, Lee SI, Lee Y, Kim EC (2009) Effect of radiopaque Portland cement on mineralization in human dental pulp cells. *Oral Surgery Oral Medicine Oral Pathology Oral Radiology and Endodontology* 108, e82–6.
- Modena KC, Casas-Apayco LC, Atta MT et al. (2009) Cytotoxicity and biocompatibility of direct and indirect pulp capping materials. *Journal of Applied Oral Science* 17, 544–54.
- Moghaddame-Jafari S, Mantellini MG, Botero TM, McDonald NJ, No' JE (2005) Effect of ProRoot MTA on pulp cell apoptosis and proliferation in vitro. *Journal of Endodontics* 31, 387–91.
- Mohammadi Z, Dummer PMH (2011) Properties and applications of calcium hydroxide in endodontics and dental traumatology. *International Endodontic Journal* 40, 697–730.
- Nekoofar MH, Adusei G, Sheykhezade MS, Hayes SJ, Bryant ST, Dummer PMH (2007) The effect of condensation pressure on selected physical properties of mineral trioxide aggregate. *International Endodontic Journal* 40, 459–61.
- Okiji T, Yoshida K (2009) Reparative dentinogenesis induced by mineral trioxide aggregate: a review from the biological and physicochemical points of view. *International Journal of Dentistry* 2009, 1–12.
- Parirokh M, Torabinejad M (2010a) Mineral trioxide aggregate: a comprehensive literature review-Part I: chemical, physical, and antibacterial properties. *Journal of Endodontics* 36, 16–27.
- Parirokh M, Torabinejad M (2010b) Mineral trioxide aggregate: a comprehensive literature review-Part III: clinical applications, drawbacks, and mechanism of action. *Journal of Endodontics* 36, 400–13.
- Pekkan G, Saridag S, Beriat NC (2011) Evaluation of the radiopacity of some luting, lining and filling dental cements. *Clinical Dentistry and Research* 35, 2–9.
- Prati C, Parrilli AP, Fini M, Dummer PMH, Gandolfi MG (2011) Interface porosity of human root canals sealed with a new flowable MTA-based sealer: a high-resolution microcomputed tomography study. *International Endodontic Journal* 44, 28.
- Rashid F, Shiba H, Mizuno N et al. (2003) The effect of extracellular calcium ion on gene expression of bone-related proteins in human pulp cells. *Journal of Endodontics* 29, 104–7.
- Sawicki L, Pamejjer CH, Emerich K, Adamowicz-Klepalska B (2008) Histological evaluation of mineral trioxide aggregate and calcium hydroxide in direct pulp capping of human immature permanent teeth. *American Journal of Dentistry* 1, 45–50.
- Schroder U (1985) Effects of calcium hydroxide-containing pulp-capping agents on pulp cell migration, proliferation, and differentiation. *Journal of Dental Research* 64, 541–8.
- Shen Q, Sun J, Wu J, Liu C, Chen F (2010) An in vitro investigation of the mechanical-chemical and biological properties of calcium phosphate/calcium silicate/bismutite cement for dental pulp capping. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 94, 141–8.
- Shubich I, Miklos FL, Rapp R, Draus FJ (1978) Release of calcium ions from pulp-capping materials. *Journal of Endodontics* 4, 242–4.
- Suh B, Cannon M, Yin R, Martin D (2008) Polymerizable dental pulp healing, capping, and lining material and method for use. International Patent A61K33/42: A61K33/42 Application number WO2008US54387 20080220; Publication number WO2008103712 (A2); Publication date 2008-08-28.
- Taddei P, Tinti A, Gandolfi MG, Rossi PL, Prati C (2009) Ageing of calcium silicate cements for endodontic use in simulated body fluids: a micro-Raman study. *Journal of Raman Spectroscopy* 40, 1858–66.
- Takita T, Hayashi M, Takeichi O et al. (2006) Effect of mineral trioxide aggregate on proliferation of cultured human dental pulp cells. *International Endodontic Journal* 39, 415–22.
- Tamburic SD, Vuleta GM, Ognjanovic JM (1993) In vitro release of calcium and hydroxyl ions from two types of calcium hydroxide preparation. *International Endodontic Journal* 26, 125–30.
- Tay FR, Pashley DH (2008) Guided tissue remineralization of partially demineralised human dentine. *Biomaterials* 29, 1127–37.
- Torabinejad M, Parirokh M (2010) Mineral trioxide aggregate: a comprehensive literature review 2013 part ii: leakage and biocompatibility investigations. *Journal of Endodontics* 36, 190–202.
- Torabinejad M, White DJ (1995) US Patent Number 5,769,638.
- Torabinejad CD, Moe H, Howley TP (1983) The effect of calcium hydroxide on porcine pulp fibroblasts in vitro. *Journal of Endodontics* 9, 131–6.
- Tronstad L, Andreasson JO, Hasselgren G, Kristerson L, Riis I (1981) pH changes in dental tissues after root canal filling with calcium hydroxide. *Journal of Endodontics* 7, 17–21.
- Tuna D, Olmez A (2008) Clinical long-term evaluation of MTA as a direct pulp capping material in primary teeth. *International Endodontic Journal* 41, 273–8.
- Wilson AD (1976) Specification test for the solubility and disintegration of dental cements: a critical evaluation of its meaning. *Journal of Dental Research* 55, 721–9.

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The Journal of Dental Research & Applications, 2007; Vol. 1(1): pg. 58-63.

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1. Ranly DM. Pulpotomy therapy in primary teeth: new modalities for old relationales. *Pediatr Dent*. 1994;16:403-409.
2. Torabinejad M, Watson TF, Pitt Ford TR. Sealing ability of a mineral trioxide aggregate when used as a root end filling material. *J Endod*. 1993;12:591-595.
3. Torabinejad M, et al. Physical and chemical properties of a new root-end filling material. *J Endod*. 1995;21(7):349-353.
4. Lee S, Monsef M, Torabinejad M. Sealing ability of a mineral trioxide aggregate for repair of lateral root perforations. *J Endod*. 1993;19(11).
5. Aqrabawi J. Sealing ability of amalgam, Super EBA cement, and MTA when used as retrograde filling materials. *Br Dent J*. 2000;188(5):266-268.
6. Torabinejad M, Ford TR, Abedi HR, et al. Tissue reaction to implanted Super EBA and mineral trioxide aggregate in the mandible of guinea pigs: a preliminary report. *J Endod*. 1995;21(11):569-571.
7. Koh ET, Torabinejad M, Pitt Ford TR, et al. Mineral trioxide aggregate stimulates a biological response in human osteoblasts. *J Biomed Mater Res*. 1997;37(3):432-439.
8. Moretton TR, Brown CE Jr, Legan JJ, Kafrawy AH. Tissue reactions after subcutaneous and intrasosseous implantation of mineral trioxide aggregate and ethoxybenzoic acid cement. *J Biomed Mater Res*. 2000;52(3):528-533.

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9. Kettering JD, Torabinejad M. Investigation of mutagenicity of mineral trioxide aggregate and other commonly used root-end filling materials. *J Endod*. 1995;21(11).
10. Keiser K, Johnson C, Tipton DA. Cytotoxicity of mineral trioxide aggregate using human periodontal ligament fibroblasts. *J Endod*. 2000;26(5):288-291.
11. Koh ET, McDonald F, Pitt Ford TR, Torabinejad M. Cellular response to mineral trioxide aggregate. *J Endod*. 1998;24(8):543-547.
12. Schwartz RS, Mauger M, Clement DJ, Walker WA 3rd. Mineral trioxide aggregate: a new material for endodontics. *J Am Dent Assoc*. 1999;30:969-975.
13. Abedi HR et al. The use of mineral trioxide aggregate cement as a direct pulp capping agent. *J Endod*. 1996;22(4):199.
14. Holland R, Souza V, Nery MJ, et al. Reaction of rat connective tissue to implanted dentin tube filled with mineral trioxide aggregate, portland cement or calcium hydroxide. *Braz Dent J*. 2001;112:3-8.
15. Holland R, Souza V, Nery MJ, et al. Reaction of rat connective tissue to implanted dentin tubes filled with mineral trioxide aggregate or calcium hydroxide. *J Endod*. 1999;25(3):161-166.
16. Pitt Ford TR, Torabinejad M, Abedi HR, et al. Using mineral trioxide aggregate as a pulp capping material. *J Am Dent Assoc*. 1996;127(10):1491-1494.
17. Pitt Ford TR, Roberts GJ. Immediate and delayed direct pulp capping with the use of a new visible light-cured calcium hydroxide preparation. *Oral Surg Oral Med Oral Pathol*. 1991;71(3):338-342.
18. Holland R, de Souza V, Murata SS, et al. Healing process of dog dental pulps after pulpotomy and pulp covering with mineral trioxide aggregate or Portland cement. *Braz Dent J*. 2001;12(2):109-113.
19. Estrela C, Bammann LL, Estrela CR, et al. Antimicrobial and chemical study of MTA, Portland cement, calcium hydroxide paste, Sealapex and Dycal. *Braz Dent J*. 2000;11(1):3-9.
20. Saidon J, He J, Zhu Q, et al. Cell and tissue reactions to mineral trioxide aggregate and Portland cement. *Oral Surg Oral Med Oral Pathol Radiol Endod*. 2003;95(4):483-489.
21. Menezes R, Bramante CM, Letra A, et al. Histologic evaluation of pulpotomies in dog using two types of mineral trioxide aggregate and regular and white Portland cements as wound dressings. *Oral Surg Oral Med Oral Pathol Radiol Endod*. 2004;98(3):376-379.
22. Funteas UR, Wallace JA, Fochtman EW. A comparative analysis of Mineral Trioxide Aggregate and Portland cement. *Aust Endod J*. 2003;29(1):43-44.
23. De Deus G, et al. Cytotoxicity of MTA and Portland cement on human ECV 304 endothelial cells. *Int Endod J*. 2005;38(9):604.

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